

**Anti-GLP-1(7-37) and GLP-1(7-36)amide (N-terminal specific)  
 Mouse monoclonal antibody**

Subclass: IgG2a/k

CAT. NO.

**ABS 033-10**

Clone:10

**SPECIFICITY** ABS 033-10 binds the free N-terminus of GLP-1(7-37) and GLP-1(7-36)amide and shows <0.2% cross-reactivity with GLP-1(9-36)amide, glucagon and human GIP. Cross-reacts 0.25% with GLP-1(1-37), approximately 6% with human GLP-2 and approximately 50% with exendin-4. ABS 033-10 cross-reacts with Liraglutide.

**IMMUNOGEN** Synthetic GLP-1(7-17)

**TESTED APPLICATIONS** ELISA, WB

**SPECIES REACTIVITY (POSITIVE)** Human

**SPECIES REACTIVITY (NEGATIVE)** Not determined

**EPITOPE SPECIFICITY** Free N-terminus of GLP-1(7-37) and GLP-1(7-36)amide.

**PRESENTATION**

**Content:** Available in 400 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.

**Preparation:** Protein-A purified

**Form:** Liquid

**Solvent:** 0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide

**Storage:** 4-8°C without exposure to light. No precautions necessary during handling.

**APPLICATION**

**ELISA:** ABS 033-10 binds free GLP-1(7-37) and free GLP-1(7-36)amide in solution. ABS 033-10 reacts in ELISA with GLP-1(7-37) coated directly onto the microtiter well.

ABS 033-10 binds free exendin-4 in solution and can function as a capture antibody for exendin-4 in combination with ABS 012-35B in order to measure Exendin-4.

When measuring "active" GLP-1 ABS 033-10B works as a detection antibody in combination with HYB 147-12.

ABS 033-10B can be used as a detection antibody in combination with ABS 046-03 as a capture antibody in order to measure non-amidated GLP-1 forms and cross-reacting about 5% with C-terminally amidated GLP-1. Results show detection limits of 44pmol/L which is 10-20 times higher than the basal concentration of GLP-1, so the assays have to be optimized.

**WB:** ABS 033-10 can be used in Western blotting.

**TARGET**

Glucagon-like peptide 1(7-36)amide (GLP-1(7-36)amide) is the principal active form of GLP-1, the other being GLP-1(7-37). GLP-1 is a peptide hormone of the glucagon family, produced by the L cells of the intestinal mucosa from the same prohormone as glucagon. The active forms are potent stimulators of glucose-dependent insulin secretion. The sequence of GLP-1 is fully conserved in all mammalian species examined so far.

**CONDITIONS**

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.