

Anti-M-Ficolin (human)**Mouse monoclonal antibody**

Subclass: IgG1/k

CAT. NO.

ABS 036-01

Clone: 1

SPECIFICITY	ABS 036-01 reacts with M-Ficolin.
IMMUNOGEN	Recombinant human M-Ficolin (aa 95 to 314)
TESTED APPLICATIONS	ELISA, IHC
SPECIES REACTIVITY (POSITIVE)	Human
SPECIES REACTIVITY (NEGATIVE)	Not determined
EPITOPE SPECIFICITY	Not determined

PRESENTATION

Content:	Available in 400 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.
Preparation:	Protein-A purified
Form:	Liquid
Solvent:	0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide
Storage:	4-8°C without exposure to light. No precautions necessary during handling.

APPLICATION **ELISA:** In ELISA ABS 036-01 reacts with recombinant M-Ficolin coated directly onto the microtiter well.
IHC: ABS 036-01 can be used in IHC. Please consult www.proteinatlas.org

TARGET M-Ficolin (Ficolin-1) is a non-serum protein that is expressed in leukocytes and lung. M-Ficolin precipitates with mannose-binding lectin (MBL)-associated serine proteases (MASP)-1 and MASP-2, indicating that M-Ficolin forms complexes with MASP-1 and MASP-2, recognizes certain types of bacteria and subsequently activates the lectin pathway.

REFERENCES**CONDITIONS**

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.