

**Anti-Oxyntomodulin (human)**

Subclass:

**Mouse monoclonal antibody**

Clone: 23

CAT. NO.

**ABS 064-23**

---

SPECIFICITY	ABS 064-23 is specific for C-terminus of human oxyntomodulin. Does not cross-react with glucagon.
IMMUNOGEN	Synthetic oxyntomodulin 31-37 coupled to carrier.
TESTED APPLICATIONS	ELISA
SPECIES REACTIVITY (POSITIVE)	Human, mouse and rat.
SPECIES REACTIVITY (NEGATIVE)	Not determined.
EPITOPE SPECIFICITY	ABS 064-23 has a different epitope specificity compared to ABS 059-20 and ABS 064-46.

---

**PRESENTATION**

Content:	Available in 400 µL and 1 ml size. 1 mg/ml +/- 15%. See Certificate of Analysis for details.
Preparation:	Protein-A purified
Form:	Liquid
Solvent:	0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide
Storage:	2-8°C

---

**APPLICATION**

**ELISA:** ABS 064-23 has been shown to work very well on a direct coat in ELISA.

**TARGET**

Oxyntomodulin is a peptide hormone of 37 amino-acid residues (29-amino acid sequence of glucagon, with an 8-amino acid carboxy-terminal extension) produced by the L cells of the gut. It modulates gastric acid secretion from the oxyntic cells of the stomach, inhibits appetite and has glucose lowering effects. Oxyntomodulin is produced by cell-type-specific processing of proglucagon that forms as well glicentin, GLP-1 and GLP-2 in the L cells of the intestinal mucosa.

**REFERENCES**

**CONDITIONS**

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.