

Anti-Vitronectin (human, dog, cat, goat, horse, bovine)

Mouse monoclonal antibody

Subclass: IgG1/k

CAT. NO.

CSI 003-23

Clone: HV23

SPECIFICITY	CSI 003-23 is highly specific for vitronectin. There is no evidence for cross-reactivity with other connective tissue proteins (fibronectin, elastin, collagen, laminin). CSI 003-23 binds nearly as well to native vitronectin as to denatured. CSI 003-23 is a potent inhibitor of integrin-mediated cell adhesion to vitronectin and a moderate inhibitor of PAI-1 binding.
IMMUNOGEN	Human vitronectin purified from plasma by heparin-affinity chromatography
TESTED APPLICATIONS	ELISA, WB
SPECIES REACTIVITY (POSITIVE)	Human, dog, cat, goat, bovine, horse
SPECIES REACTIVITY (NEGATIVE)	Not determined
EPITOPE SPECIFICITY	Epitope is located in the somatomedin B domain

PRESENTATION

Content:	Available in 400 µL and 1 mL size. 1 mg/mL +/- 15%. See Certificate of Analysis for details.
Preparation:	Protein-A purified
Form:	Liquid
Solvent:	0.01 M phosphate buffer, pH 7.4, containing 0.5 M NaCl and 15 mM sodium azide
Storage:	4-8°C without exposure to light. No precautions necessary during handling.

APPLICATION	ELISA: CSI 003-23 also binds to vitronectin in ELISA when vitronectin is coated directly onto the microtiter well. In Western blotting a dilution guideline of 1/50 and 1/200 has proved successful (1, 2). WB: CSI 003-23 can be used in Western blotting (1).
-------------	--

TARGET	Vitronectin is a plasma glycoprotein that circulates in the blood. Vitronectin is circulating as a mixture of both 75 kDa and 65 kDa forms. Vitronectin is a major cell adhesive glycoprotein and is a common component of extracellular matrix and plasma. It competes effectively with other plasma proteins and is often involved in cell attachment, regulation of blood coagulation and immune responses. It has similar tissue distribution to fibronectin and also its integrin receptor recognizes fibronectin (2).
--------	---

REFERENCES	<ol style="list-style-type: none">Morris CA, Underwood PA, Bean PA, Sheehan M, Charlesworth JA (1994) Relative topography of biologically active domains of human vitronectin. Evidence from monoclonal antibody epitope and denaturation studies. J Biol Chem 269:23845-23852.Underwood PA, Kirkpatrick A, Mitchell SM (2002) New insights into heparin binding to vitronectin: studies with monoclonal antibodies. Biochem J 365:57-67.
------------	--

CONDITIONS

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.