

Anti-Placental protein 14 (human, PP14, glycodelin A)

CAT. NO. BTE 001-16

OVERVIEW

Product Name	Anti-Placental protein 14 (human, PP14, glycodelin A)	Conjugation	Unconjugated
Description	Mouse monoclonal antibody	Host	Mouse
Isotype	IgG1/k	Clone	16
Tested Applications	ELISA, IHC		

SPECIFICITY

Specificity	BTE 001-16 is specific for human PP14 protein 14 (human, PP14, glycodelin A). Epitope differs from BTE 001-13		
Immunogen	Placental protein 14 (glycodelin A) purified from second trimester amniotic fluid	Gene ID	5047
Target	Human placental protein 14 (PP14; also known as glycodelin and progesterone-associated endometrial protein) is a protein of the lipocalin structural superfamily. PP14 is the most abundant product of the secretory endometrium, and has been proposed as a biochemical marker of endometrial function in women undergoing in vitro fertilization treatment. (1)		
Species Reactivity POSITIVE	Human	Species Reactivity NEGATIVE	Not determined

PROPERTIES

Form	Liquid	Unit Size	0,4 mL and 1 mL
Concentration	1 mg/mL \pm 15%, See CoA for lot details		
Purification	Protein A or Protein G purified	Purification Notes	BSA free
Storage buffer	0.01 M phosphate buffer, pH 7.4, with 0.5 M NaCl and 15 mM sodium azide		
Storage condition	2-8°C without exposure to light		
Safety	Wear protective clothing		

TESTED APPLICATIONS

ELISA	BTE 001-16 reacts strongly to human PP14 in sandwich ELISA in combination with a polyclonal antibody.
IHC	When staining formalin-fixed paraffin-embedded endometrial tissue from the late secretory phase, staining is restricted to the glandular cells. (2+3)

SCIENTIFIC REFERENCES

- Westergaard LG, Wiberg N, Andersen CY, Laursen SB, Kliem A, Westergaard JG, Teisner B (1998) Circulating concentrations of placenta protein 14 during the natural menstrual cycle in women significantly reflect endometrial receptivity to implantation and pregnancy during successive assisted reproduction cycles. Hum Reprod 13:2612-2619.
- Hustin J, Philippe E, Teisner B, Grudzinskas JG (1994) Immunohistochemical localization of two endometrial proteins in the early days of human pregnancy. Placenta 15:701-708.
- Tornehave D, Fay TN, Teisner B, Chemnitz J, Westergaard JG, Grudzinskas JG (1989) Two fetal antigens (FA-1 and FA-2) and endometrial proteins (PP12 and PP14) isolated from amniotic fluid: localisation in the fetus and adult female genital tract. Eur J Obstet Gynecol Reprod Biol 30:221-232.

CONDITIONS

Unless otherwise marked, all products are for research use only. Not for use in diagnostic procedures. Not for use in human therapeutic applications. For in vitro use or further manufacture only. The information and product are offered without guarantee as the ultimate conditions of use are beyond our control. The foregoing is in lieu of all warranties, expressed or implied, including implied warranties of merchantability and fitness for a particular purpose. In no event shall BioPorto Diagnostics A/S be responsible for loss of profits or indirect consequential losses resulting from use of its products.